



# MALDI Oligonucleotides Sample Preparation Protocol

---

## 1. Plate Type



MALDI plate Ground Steel



MALDI plate Small Anchorchip (800µm)

## 2. Plate Cleaning

### Cleaning solvents:

Solvent A: 2-Propanol

Solvent B: Water

Solvent C: TA30

### Cleaning protocol:

- Wipe off old preparations using tissue wetted with solvent A and solvent B, subsequently.
- Sonicate the target plate for 10min using solvent A.
- Sonicate the target plate for another 10min using solvent C.
- Dry the target plate using high purity nitrogen or air. Don't wipe the target front side of the cleaned plate. If pure gases are not available, let dry the plate under room conditions.

## 3. Matrix solution preparation according to plate type

### 3.1 Plate type: Ground steel

- **Matrix Solution:** 3-HPA saturated in H<sub>2</sub>O:ACN 50:50 (v:v), containing 10mg/ml Diammonium hydrogen citrate ( DAC)

### 3.2 Plate type: AnchorChip (800µm)

- **Matrix Solution:** 3-HPA 10mg/ml in H<sub>2</sub>O:ACN 50:50 (v:v), containing 10mg/ml DAC

### 3. Sample spotting

- **Sample Solvent:** water
- **Preparation protocol:** Apply 0.5 $\mu$ l of the matrix solution on the target. Let the matrix dry down at room temperature. Add 0.5 $\mu$ l of the sample solution on top of the dried matrix preparation spot. Let the preparation dry down at room temperature.